

MEDISTIM VERIQ™ SYSTEM USER MANUAL

Mod. VQ4122 | VQ2111 | VQ1111 | VQ4001 | VQ2011 | VQ1011 | VQ1001 |



CONTENTS

1. INTRODUCTION	4	5.4 Preparation by surgeon	12
1.1 Purpose	4	5.5 Tasks during procedure by operator	13
1.2 Contact information	4	5.5.1 Doppler measurements	13
1.3 User Manual structure and use	4	5.5.2 Performing vessel search	13
1.3.1 Relevant books and documents	4	5.5.3 Detect position and quantify a stenosis	13
1.3.2 User Manual Update	4	5.5.4 Transit time flow measurements	13
1.4 Warning statements and safety markings	4	5.5.5 Optimizing the Doppler and flow display	13
1.4.1 User Manual warning statements	4	5.5.6 Saving Measurements	14
1.4.2 User Manual language options	4	5.5.7 Memory scroll and trace review	14
1.5 Safety marking description	4	5.5.8 Review of recorded data	14
2. INDICATIONS FOR USE	5	5.6 Tasks during operation by surgeon	15
2.1 Indications for use statement	5	5.6.1 Doppler velocity measurements	15
2.2 Contraindications	5	5.6.2 Flow measurements	15
2.3 Health, environment and safety	5	5.7 Post-operative procedures	16
2.4 Warranty	5	5.7.1 Retrieve patient from database	16
2.5 Operational environment requirements	5	5.7.2 Cleaning and disinfection of probes	16
2.6 Required qualifications of personnel	5	5.7.3 Sterilization	16
2.7 Limitations for use	5	5.7.4 Reporting	16
2.8 Possible operational influence	6	5.8 Icon description	17
2.9 Consumables	6	5.9 System settings description	18
2.10 Access to spare parts and service	6	5.10 Other adjustments	18
3. TECHNICAL DESCRIPTION	7	5.10.1 Tuning Doppler noise	19
3.1 Front panel and channel configurations	7	5.11 Data Management	20
3.1.1 Pay-Per-Procedure	7	5.11.1 Search function	20
3.1.2 USB port connections	7	5.11.2 Edit, print, delete and export measurements	20
3.2 Principle of measurements	7	5.11.3 Export and import data	20
4. DELIVERY AND CONFIGURATION	8	5.11.4 Safe removal of storage devices	20
4.1 Receiving the Medistim VeriQ™ system	8	5.12 Derived traces	21
4.2 Unpacking and packing	8	5.12.1 Fast Fourier Transformation (FFT)	21
4.3 Verifying system identity	8	5.13 Trend measurements	22
4.4 Verifying electrical safety	8	5.13.1 Save, review and edit trend curves	22
4.4.1 Starting the system	8	5.14 Measurement and Analysis (M&A)	22
4.5 Initial configuration	8	6. TROUBLESHOOTING	24
4.6 Auxiliary inputs	9	6.1 Service of the VeriQ™ System	24
4.6.1 Calibrating the AUX channel for ECG	9	6.2 Basic troubleshooting	24
4.6.2 Calibrating the AUX channel for pressure	9	6.2.1 Troubleshooting users guide	24
4.6.3 Shutting down the system	9	7. MAINTENANCE	25
5. OPERATING THE VeriQ™	10	7.1 Probe statistics	25
5.1 Requirements for personnel	10	7.2 Periodic maintenance	25
5.2 Basic principles of operation	10	7.3 Storage	25
5.2.1 Correct plugging and unplugging of probes	10	7.4 Packing and transport	25
5.3 Preparations by system operator	11	7.5 Disposal instructions	25
5.3.1 Probe preparation	11	8. TECHNICAL SPECIFICATIONS	26
5.3.2 System Start-up	11	8.1 VeriQ™ System overview	26
5.3.3 Entering patient data	11	8.2 Technical specifications	27
5.3.4 Connecting ECG (optional)	11	APPENDIX A. PROBE APPLICATION OVERVIEW	31
5.3.5 Inserting key-card (PPP systems only)	11	A.1 Probe application and size guide	31
5.3.6 Probe verification	11	A.2 List of available Medistim probes	32
5.3.7 Defining the vessel names	12		

1. INTRODUCTION

1.1 Purpose

The purpose of this User Manual is to give a thorough description of the Medistim VeriQ™ System.

1.2 Contact information

Medistim ASA
Head Office:
Økernveien 94
0579 Oslo, Norway

Telephone: +47 23 05 96 60
Fax: +47 23 05 96 61
E-mail: medistim@medistim.com
Internet: <http://www.Medistim.com>

Manufacturing and Service support:
Moloveien 10,
3187 Horten, Norway

Telephone: +47 33 03 17 10
Fax: +47 33 03 17 26
E-mail: service@Medistim.com

1.3 User Manual structure and use

1.3.1 Relevant books and documents

Please refer to the following relevant information in addition to this User Manual:

- www.medistim.com
- Service Manual for VeriQ™p/n: VQ1980GB
- "Intraoperative Graft Patency Verification during On- and Off-Pump Coronary Bypass Surgery" by Hratch L. Karamanoukian, MD and Harry W. Donias, M.D., Medistim ref. number MM117803
- "Intraoperative Graft Patency Verification in Cardiac and Vascular Surgery". ISBN# 0-87993-488-3s

1.3.2 User Manual Update

In the incidence of additional modifications to the VeriQ™, a new revision of the Medistim VeriQ™ User Manual will be distributed by Medistim ASA.

1.4 Warning statements and safety markings

1.4.1 User Manual warning statements

This User Manual uses Warnings and Notes according to the following definitions:

Warning

This warning will describe clinical contraindications and possible damage to the device if the recommended instructions or recommendations are not followed. Please read and follow

these warnings carefully. If Warnings are not taken into consideration, serious injury and damage to personnel, devices and the operating environment can occur. The manufacturer cannot be held liable for injury or damage if these precautions are disregarded.

Note

A Note contains important tips, recommendations and supplementary information intended to optimize the use of the system.

1.4.2 User Manual language options

The User Manual is, upon request, available in English, French, Spanish, Swedish, Danish, Norwegian, German and Italian.

1.5 Safety marking description



On / Off button (push-push)



Attention, consult accompanying documents



Protective earth (ground)



Terminal for Potential Equalization



Consult accompanying documents
(Required to consult for safety)



Pushing prohibited

Warning

Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the instructions provided in this manual. Portable and mobile radio frequency (rf) communications equipment can affect medical electrical equipment. No connection from the internal signal input/outputs located at the lower side of the electronic unit to external equipment are allowed. For DICOM enabled systems, the network connection marked with IEC60601-1 compliance are excepted from this.

Warning

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth

Warning

Equipment connected to the systems signal input/ output shall comply with IEC 60601-1 clause 16

Potential Equalization conductor

A terminal for potential equalization is located at the back of the power unit of the system. The terminal can be used if potential equalization needed but is no replacement for the protective earth conductor of the power cord as required by IEC60601-1 cl. 8.6.7.

2. INDICATIONS FOR USE

The Indications for use statement for the Medistim VeriQ™ System is shown in Section 2.1 below. For the System model that supports Ultrasound Imaging, VeriQ C™, refer to indications for use statement in: “Medistim VeriQ C™ System User Manual, Addendum to the Medistim VeriQ™ System User Manual VQ1990 for model VQ4122C.”

2.1 Indications for use statement

The Medistim VeriQ™ System is designed to perform intra-operative guidance and quality control during cardiovascular surgical procedures and meets the demands for documentation of surgical procedures. The following is the indication for use statement cleared by the US Food and Drug Administration on April 29th 2004 (K040228):

“The Medistim VeriQ™ System is an intraoperative diagnostic system that utilizes ultrasonography to guide surgeons to successfully plan and accomplish surgical interventions. The clinical indications for the device are:

1. Accurate transit time blood volume and Doppler velocity flow measurements during cardiovascular-, vascular-, transplantation- and neurosurgery.
2. Simultaneous measurements of blood pressure, vascular resistance, interfaced physiological signals and other derived parameters during these procedures.
3. Detection of normal and abnormal blood volume and Doppler velocity flow patterns during these procedures.
4. Provides guidance to prepare surgical plans at the initiation of surgery and to support the successful accomplishment of surgery including detection and location of vessels during surgical procedures.
5. Detection and quantification of the degree of stenosis in arteries by using the Doppler velocity profile.”

2.2 Contraindications

Warning

Do not use this device for other applications than the intended use. The Medistim VeriQ™ System is not intended for any kind of fetal applications.

Warning

Transit time and Doppler flow measurements are routinely used for non-diseased vessels such as saphenous vein grafts, different arteries and harvested mammaries. In the vessels, thrombus or calcified plaque will not be present. However, thrombus or calcified plaque may be present before certain clinical procedures, such as carotid endarterectomy, or before performing a proximal anastomosis on the aortic arch. In such, or similar procedures the surgeon may want to measure blood volume flow or velocities prior to performing the surgical intervention.

Manipulation of vessels with thrombus, calcification or other conditions may have serious consequences for the patient, and represents a contraindication for the Measurement of

blood flow. Sound medical judgment should therefore be exercised and the final responsibility for measuring flows in diseased vessels lies with the physician. Please avoid any pinching of vessels during the flow measurements.

2.3 Health, environment and safety

No health, environment or safety considerations are related to the Medistim VeriQ™ System.

2.4 Warranty

The Medistim VeriQ™ system is warranted against defective material and poor workmanship for a period of 24 months after installation, and limited to 25 months from date of shipment.

The Medistim probes are warranted against defective material and poor workmanship for a period of 12 months after first use, and limited to a maximum of 50. (See specifications for each probe.)

All other probes are warranted against defective material and poor workmanship for a period of 6 months after first use, and limited to a maximum of 50 usages.

2.5 Operational environment requirements

The following are requirements for use:

- Correct electrical supply
- Ground terminal
- Indoor use only
- For intra-operative use only

2.6 Required qualifications of personnel

All persons operating VeriQ™ must have read this User Manual before using the device.

2.7 Limitations for use

Transit Time Flow Measurement (TTFM) and Doppler flow measurements are routinely used for non-diseased vessels such as saphenous vein grafts, different arteries and harvested mammaries. In these vessels, thrombus or calcified plaque will not be present.

However, thrombus or calcified plaque may be present before certain clinical procedures such as carotid endarterectomy or before performing a proximal anastomosis on the aortic arch. In such procedures, or similar procedures, the surgeon may want to measure blood volume flow or velocities prior to any surgical intervention.

Manipulation of vessels with thrombus, calcification or other malfunctions in order to measure blood flow may have serious consequences for the patient and is a possible contraindication for performing blood flow measurements. Sound medical judgment should therefore be exercised and the final responsibility for measuring flows in diseased vessels will lie with the physician. Any pinching of the vessels during flow measurement should be avoided.

2.8 Possible operational influence

As probe size, positioning and placement on the vessel may influence accuracy, it is important that all measurements are made as described in this manual. To avoid the disruption or corruption of the calculated Diastolic Filling Percentage (DF), no action should be taken that compromises the quality of the ECG signal, and DF is dependent on a stable ECG recording. Simultaneous use of Diathermy and flow/velocity measurements should be avoided as the Diathermy interferes with the measurements and can make them unreliable.

2.9 Consumables

Necessary consumable items to perform an operation:

- Various probes
- Sterile Ultrasound couplant
- Ink cartridges for color printer
- Printer paper

2.10 Access to spare parts and service

Please contact your local representative to request VeriQ™ service or spare parts.

Refer to our website www.medistim.com for a list of local contacts and distributors.

Warning

No modification of this equipment is allowed.

Warning

Do not modify this equipment without authorization of the manufacturer.

Warning

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

3. TECHNICAL DESCRIPTION

The VeriQ™ system is designed for intra-operative quality control and the need for documentation of surgical procedures. Through its capability to measure transit time volume flow and Doppler flow velocity, the VeriQ™ can help develop optimal grafting strategies and to document graft patency.

The VeriQ™ system consists of a 19 inch touch screen, displaying significant flow variables. The menu-driven program is managed through the use of an intuitive graphical, icon-driven interface. Improper data entry and incorrect operations are minimized through the use of an on-screen keyboard that appears only when it is required.

The VeriQ™ uses five main screen views; the Home screen, the Measure screen, the Data screen, the Help screen and the System Settings screen. The user can move between screens by pressing the screen tabs located in the top left corner. All reports and documentation can be written directly to the printer connected to the system.

3.1 Front panel and channel configurations

The VeriQ™ system can be supplied with a range of different channel configurations to accommodate various clinical needs. The channel configurations are described in the model number, where the first digit lists the number of flow channels, the second digit is the Doppler channel, third digit is the number of pressure channels, and the fourth digit is the number of auxiliary channels.



The VeriQ™ 4122 has four flow channels, one Doppler channel, two pressure channels and two auxiliary channels. The model number is part of the system identity and can be found on the front of the system under the blue top cover.

3.1.1 Pay-Per-Procedure

The VeriQ™ system is either supplied as a capital purchase, or through a pay-per-procedure (PPP) program, where license “smart cards” are required in order to activate the device.

3.1.2 USB port connections

The VeriQ™ system has 2 USB 2.0 ports located on the left side of the front panel for connecting external equipment like a memory stick to import or export data.

Warning


The front panel USB connector has no galvanic separation. In order to maintain continued patient safety, no equipment powered by an external power source must be connected to the

USB port during clinical use of the system.

3.2 Principle of measurements

Transit time flow volume

Medistim transit time probes surround the vessel, generating a uniform ultrasound field across the vessel lumen. Ultrasound pulses are transmitted from a crystal on one side of the probe, to the reflector, and back to a crystal on the opposite side. The transit time is measured for each pulse, with the difference in transit time between the pulses going upstream and downstream being proportional to the volume (Q) passing through the probe.




Doppler principle

When an ultrasound beam is reflected by a moving object, the frequency of the reflected pulse is changed. An object moving towards the ultrasound beam will compress the waveform and increase the frequency. Correspondingly, an object moving away from the beam will lengthen the waveform and decrease the frequency.

The change in frequency, also called Doppler shift, represents the velocity and direction of the moving target. To limit problems relating to the ultrasound beam's angle of incidence with the blood, the Doppler probe has a built-in angle of 45 degrees. Holding the probe perpendicular to the vessel direction will ensure accurate velocity measurements. Note that this causes the probe to measure the velocities in front or behind the actual probe position, but never right under the probe.

The VeriQ™ applies a pulsed Doppler, allowing the user to control the depth from where the velocity should be measured. Additionally, sample volume is adjustable, and selected as a range around the depth setting. The default settings are a depth of 5 mm, and a volume of 6 mm. These settings allow for sampling flow velocities at depths from 2 to 8 mm from the probe surface.



4. DELIVERY AND CONFIGURATION

The purpose of this section is to prevent damage to property and injury to personnel through the delivery, installation and preparation of the VeriQ™.

4.1 Receiving the Medistim VeriQ™ system

The system is shipped complete in one large, wooden crate. Probes and accessories may be packed in a carton placed inside the system crate.

Upon receipt, please observe that there is no visible damage to the crate, and that the handling indicators (shock and tilt) attached to the crate are still intact.

Any visible damage should be reported to the freight forwarder upon receipt.

4.2 Unpacking and packing

Please refer to the instructions provided with each crate for proper unpacking.

Medistim recommends that the original packing materials be used when the system is transported between facilities. For service purposes, a specially-designed container is available through your local representative.

4.3 Verifying system identity

The model number, serial number and revision number are noted on the system identity label located on the upper front side of the electronic unit, just under the blue top cover.

4.4 Verifying electrical safety

The system has been tested for electrical safety prior to delivery, and a report is attached to the system configuration sheet placed together with the User Manual in the shelf under the printer.

If a new verification is made in connection with the installation, please refer to the chapter 8.2, External wiring, for system grounding information.

Warning

The Medistim VeriQ™ System can be seriously damaged if the system is turned off before the Windows booting sequence is complete. Do not turn the device off during the booting sequence. Wait until the patient data entry screen is displayed.

4.4.1 Starting the system

Check to make certain the main power switch, located on the right hand side underneath the printer shelf, is in the “on” position. When this switch is on, both the switch itself and the printer indicator lights will be on. To start the VeriQ™, press the soft-power button, located on the right side of the top cover beside the monitor. The software boot sequence will take approximately 35 seconds before the system is ready for use.

4.5 Initial configuration

Start the VeriQ™ and select the “system settings” screen. The system setup menu is accessed by pressing the System settings icon in the top left of the screen.

Software version

The software version currently installed on the system can be identified by opening the “System Information” page located on the “Advanced System Functionality” page.

Enter hospital name and address

The hospital name and address should be entered before the system is put into service. On the system settings screen press “edit” in the upper left, and a touch screen keyboard will appear. Enter the name and address of the hospital in the available dialog boxes. Press OK when completed.

Select system profile

The VeriQ™ can handle different profiles; Cardiac, Vascular and Experimental. adapting the available vessel names to the selected profile. Each profile can be edited, and operating surgeons can be stored under each profile. Select the appropriate profile by pressing the drop-down menu “System profile”

Select ECG Channel

The VeriQ™ should always be hooked up to the anesthesia monitoring system in order to slave the ECG signal. The ECG signal is used to synchronize the timing of the recorded flow or velocity traces. VeriQ™ will be supplied with an ECG cable, fitting one of the two AUX connectors, and should be matched to the manufacturer of the anesthesia monitoring system. Verify that the AUX channel selected for ECG in the system setup menu corresponds to the available ECG cable.

Set measuring unit

Metric or US measuring units are selected by pressing the drop down menu “Measurement unit”. Note that the measuring unit is for display use only, and can be changed at any time without influencing data already stored in the system.

After connecting a pressure transducer, the system allows the setting of zero-level and fine-tuning of the gain for optimal measurements. The procedure requires reference signals of 0 mmHg and 100 mmHg. Enter the calibration screen by pressing the Pressure properties button.

First apply 0 mmHg and press the zero-button. Wait until the curve stabilizes at zero. Then apply 100 mmHg and adjust the gain until the curve rests at 100 mmHg.

The calibration of the pressure channels are non-persistent and must be done for each new transducer.

4.6 Auxiliary inputs

The auxiliary inputs can be used to measure any external voltage signal within $\pm 4V$. The signal is filtered, conditioned, assigned to a measurement unit and included in the VeriQ™ measurement set. The available signal ranges are $\pm 10mV$, $\pm 100mV$, $\pm 500mV$ and $\pm 4V$. The user can select signal level, measurement unit, filtering and signal condition from the AUX calibration screen as described in the sections below. When selecting mmHg as the measurement unit, the system treats the AUX input as a pressure input with respect to derived traces and calculations.

The AUX channel selected for ECG input will automatically display the heart rate in BPM and the systolic and diastolic phases can be drawn on the curve trace (user selectable).

For other measurement units selected, the curve trace will simply show the measured value in terms of the selected unit.

4.6.1 Calibrating the AUX channel for ECG

Turn on the system and connect the ECG cable to the anesthesia monitoring system. For the purpose of calibration, an ECG simulator can be used with the monitoring system. Select the **Live session** screen, and the **Auxiliary curve** window will appear. Depending on signal strength, the ECG waveform will be visible. If the scale and measuring unit appear correct, no further adjustments are necessary.

To improve the settings for optimal display, press the **ECG properties** button and the **Calibrating ECG** window will appear.

After selecting the correct measuring unit, ECG, select an input factor that corresponds to the signal input. If the displayed signal has a low amplitude, then select a smaller input range. Note the AD range displayed in the upper right-hand corner of the curve window, as it is important that the AD range is less than the saturation level of 100%. For fine adjustments within the selected range, use the **Scaling up** and **down** buttons. When the correct ECG curve is displayed, the settings should be saved by pressing the **Save** button. To ensure that the new setup is used whenever the ECG is connected, press the **Set current as default** button.

Note

When using the Medistim VeriQ™ System with more than one monitoring system that requires independent adjustments, different calibration settings can be saved under a profile name, such as OT 1 and OT 2. The operator can choose the stored calibration each time the system is moved by entering the calibration menu and selecting the appropriate saved profile.

4.6.2 Calibrating the AUX channel for pressure

Turn on the system and insert the cable from the signal monitor. For calibration, a simulator that emits known signals can be used. Select the **Live session** screen, and the auxiliary curve window will appear. Depending on the signal strength, the monitored curve may be visible straight away. To calibrate the curve, press the **AUX1/2 properties** button and the **Calibrating AUX1/2** window will appear.

In order to calibrate the curve properly, a signal with two known amplitudes is needed. One of these is typically zero, the other can be chosen randomly as long as it is feasible and known.

First, select the correct measuring unit. (For pressure, it would be mmHg.) Then select an input range that corresponds to the physical signal. If the signal has low amplitude, select a smaller input range.

It is important that the AD range displayed in the upper part of the curve window is less than the saturation level of 100%. When the signal from the monitor is set to zero, press the **Zero** button on the lower right side of the curve view to set the zero signal level. Change to the signal value with known amplitude and use the **scale up** and **scale down** buttons to adjust the displayed value so that it fits the input. Readjust the input range if necessary.

When the calibration is satisfactory, assign the profile an appropriate name in the **Current profile** edit box and press **Save**.

To set the active profile as default for the channel, press the **Set current as default** button on the lower left part of the screen.

For all types of AUX measurements the calibration process is exactly the same. Exchange the unit to suit the measured variable. Adjust input range and gain so that the signal is within saturation limits and correct for two known values.

4.6.3 Shutting down the system

The system should be turned off by the soft-power button first, and be allowed to power down before using the main power switch on the transformer unit. Do not unplug the power cable before this is done.

It is not necessary to exit the application or do additional saving before powering down the system.

When preparing for storage, unplug the power cord from the wall by pulling on the plug itself, do not pull on the cable. If the system is equipped with a network cable, be careful to also unplug this from the wall socket. Use the cable hook on the back cover of the VeriQ™ to coil the cable(s). All probes should be unplugged and stored properly if not sent to cleaning and sterilization.

5. OPERATING THE VeriQ™

The purpose of this section is to provide the user with a reference for the operation of the VeriQ™.

5.1 Requirements for personnel

To successfully operate the VeriQ™, at least two persons are required. One controlling the probe and the other the system controls.

The main task for the system operator is to prepare the system before use, adjust the system during operation and to save essential findings in order to document the procedure.

During the operation, it is important for the surgeon and operator to maintain a good dialog to work effectively and get the best possible results.

5.2 Basic principles of operation

VeriQ™ incorporates 19-inch touch screen with a graphical user interface. Operation of the VeriQ™ is intended to be very intuitive, eliminating non-essential commands, and presenting an on-screen keyboard only when necessary. To further simplify operation for novice users, certain advanced functions are not readily-accessible. Advanced functions can be activated once the operator is familiar with the basic system operation.

The patient name will always be displayed at the top center of the screen. When no patient data has been entered or selected, the system will use the current date and time as the patient name. In the top right corner the current data and time is displayed. On Pay-Per-Procedure (PPP) systems a field in the top right corner will display the remaining procedure time once a procedure has been activated.

The VeriQ™ has five main screen views:

1. Home screen
2. Measure screen
3. Data screen
4. Help screen
5. System Settings screen.

The user can move between screens by pressing the screen tabs located in the top left corner. The "Home" screen is used to enter patient data when starting a new surgical procedure and when searching for a patient in the database. This is the start-up page, and by pressing the "Home" tab the system will return to this screen.

The "Measure" screen is the main screen for flow measurements, and can be accessed by pressing the "Measure" screen tab or by pressing "OK" on the patient data screen. Systems installed based on a "pay-per-procedure" model must be activated by using a smart card. On-screen prompts will assist in system activation, and when a smart card has not been inserted into the system, the message "Insert valid smart card to activate" will be shown in the message bar.

System status is displayed in the message bar that appears near the top of the VeriQ™ monitor. To alert the user of status changes, the background color will change. Press the "Data" tab to enter the screen that presents a thumbnail view of each recorded curve on the current selected patient. Here, stored datasets can be selected for viewing, editing, reporting or exporting. The "Help Screen" contains an electronic User Manual, and video tutorials. The "System settings" screen contain the user configurable settings for the VeriQ™.

5.2.1 Correct plugging and unplugging of probes


To avoid accidental errors related to plugging a probe in the wrong connector, the VeriQ™ probe panel and all probes have unique keying features. This ensures that only TTFM probes can be plugged into the flow channels and vice versa for the other channels.

When plugging a probe into the system connector, hold the probe with the embossed arrows facing upwards, see picture below, and press straight inwards.



All Medistim probes require very little force to be inserted or removed. If excessive force is needed to plug a probe this is a sign that something is damaged either in the probe connector or on the system side connector. Inspect both connectors carefully to identify possible bent or broken pins in the probe connector, plugged or otherwise obstructed holes in the system side connector.

Unplugging the probe is done by gripping the probe on the widest part of the connector and pulling straight out. This will release the locking mechanism and the probe will unplug correctly. Do not attempt to unplug by pulling on the cable or strain relief, this will not release the locking mechanism and may damage the probe. See below image correct and incorrect removal technique.



5.3 Preparations by system operator

5.3.1 Probe preparation

Locate the necessary probes to perform the planned measurements. For CABG procedures, a variety of QuickFit probes will be required to verify graft patency. Additionally, a Doppler probe will be needed to detect vessels, or locate stenoses.

To select the optimal probe sizes, please refer to section 10 Probe Assortment.

5.3.2 System Start-up

Turn the system on before taking measurements, in order to interface the necessary auxiliary signals and to enter patient information.

Make sure the following items are available before the procedures:

- A selection of sterile probes suitable for the planned procedure
- The ECG interface cable
- A valid key card (PPP systems only)

5.3.3 Entering patient data

When the system is turned on, it will start in the Home screen. Use the soft-keyboard on the screen to enter patient information, select the surgeon, and press OK. The system will automatically proceed to the Measure screen.

Note

Medistim recommends that the patient name is entered in an orderly and uniform way for easy storage and retrieval, e.g. always with the last name first. If no patient name is entered, the date and time of operation will be used as the patient identifier.

5.3.4 Connecting ECG (optional)

Connect the ECG cable from the monitoring system and to the appropriate Aux connector on the VeriQ™. Verify that the ECG waveform is displayed correctly on the screen. For adjustments of the ECG signal please refer to section 4.5.3.

The system may be left in this state until the surgeon is ready to start measuring.

5.3.5 Inserting key-card (PPP systems only)

In order to do measurements, a valid license key card (Smart card) is required for PPP-enabled systems. Insert the license card into the card reader on the front of the device and await the dialog box confirming validity of the license. The procedure timer of five hours (seven hours for imaging licenses) will start when "Activate" is pressed. The ECG curve will be visible and the "Insert valid smart card to activate" message will disappear.

Medistim supplies different smart cards to activate the system. These are the different smart cards versions:

1. Single use full function (transit time and Doppler).
2. Single use flow (transit time only).
3. Single use Vascular.
4. Single use Imaging Add-On.
5. Single use flow + Imaging.
6. Time-Limited full function (transit time and Doppler).
7. Time-Limited flow + Imaging.

The "single use" licenses are disposable and require that a new card be used for every new procedure. The "time limited" licenses opens the given functionality for a period of time. At the first use, the current time is registered in the card and the countdown timer starts. During the period specified by the license there are no limitations on the number of procedures performed.

The "Imaging Add-on" license can only be used after a flow license has been activated. It activates the Imaging modality and resets the procedure timing. The procedure time for all imaging licenses is seven hours. For more information about the VeriQ™ System model(s) that supports Ultrasound Imaging, VeriQ C™ refer to: "User Manual VeriQ C™, Addendum to VeriQ™ System User Manual VQ1990 for model VQ4122C."

Note

The validity of license cards can be tested at any time by inserting them in the card reader. The system will display the status of the card, and if cancel is pressed or the card removed without pressing "Activate procedure", the card license will remain intact.

5.3.6 Probe verification

Before probes are connected to the system, the monitor will display ECG traces.


1. Connect the selected probes to the respective system channels (flow, Doppler or Imaging.)
2. For each entered probe a vessel name dialog box will appear. For the probe verification purpose the vessel name is not important, just press cancel until the measuring screen is displayed.
3. Verify that the probe identity of the connected probes appear on the screen.



Medistim recommends placing the probe head in sterile saline and connecting it to the system before the surgeon is ready to measure. This allows each probe to be tested before use and it will improve the acoustical coupling when placed on a vessel. With the probe immersed in saline, all transit time probes should obtain an ACI value of > 90%. If a lower value is obtained, shake the probe in the saline solution to remove any bubbles. Air bubbles on the probe will significantly affect the ACI value. After removing any interfering air bubbles, re-check the ACI value. If the displayed ACI value is not a minimum of 90%, the probe is not working properly and should be replaced.

Note

On PPP systems, acoustical coupling can be tested before inserting the key card. The system will only display the acoustical coupling indicator and other probe related information.



5.3.7 Defining the vessel names

When a probe is connected, a Vessel name dialog page will automatically appear.

Note

The automatic appearance of the dialog page can be switched off in the system settings menu.


During a measurement, or if the graft dialog does not automatically come up, press the button "Vessel name Undefined" to launch the graft dialog. This button will also be displaying the previous graft selection if one is made i.e. LIMA - LAD. Press this button to get into the graft dialog to set a new graft name.

Register vessel


1. Select type of graft conduit, source, target artery and any attributes in the dialog page (1). The selection can also be made from the list of the five mostly used graft combination on the current machine.(2)
2. The free text field (3) can be used when the vessel

names already given do not cover the grafts and vessels in question. Only one of these selections can be made and the present selection will exclude the previous.

The selection will be displayed under current graft directly above the OK button. To finish, press OK (4) to select or cancel to go back to the measurement without doing the change/selection.



When generating a report or reviewing the case at a later time, the vessel name will be important.



Note

The source may be omitted when obvious, and is primarily intended for sequential grafts. The top five most commonly used graft names on the current machine can be selected directly on the screen

Searching for vessels and detecting stenoses with the Doppler probe will only require a target vessel and, possibly, an attribute if measurements from different positions on the same vessel are stored.


5.4 Preparation by surgeon

The surgeon needs to judge the diameter of the vessels to be measured and inform the operator which probe(s) to connect. Arterial conduits such as the internal mammary arteries have to be skeletonized for a length corresponding to the width of the probe, approximately 1cm.

5.5 Tasks during procedure by operator

5.5.1 Doppler measurements

The VeriQ™ will display the Doppler spectrum at the default 5 seconds sweep rate as soon as a probe is connected, as shown on figure below.



Note


The Doppler velocity curve is sensitive to the flow direction. If the signal appears under the baseline (see figure above), the Doppler probe needs to be rotated 180° to accommodate the measured flow direction.

5.5.2 Performing vessel search

When searching for intramuscular vessels, the surgeon places the probe in the approximate area of the vessel, and listens for the audible Doppler signal. The operator should select an appropriate volume setting that enables the surgeon to hear the signal, and possibly help differentiate arterial flow (pulsatile waveform as the figure above) or venous flow (continuous flow) from noise originating from probe movements.

5.5.3 Detect position and quantify a stenosis

When looking for a stenosis, the surgeon first needs to measure a normal, patent segment of the vessel. When a stable curve is displayed, the operator should press the "reference button".



The system will display a reference line, demonstrating the recorded reference peak velocity. When the surgeon relocates the probe, the measured peak velocity will be compared with the reference line. The system will display the change in peak velocity as a percentage of stenosis. The velocity scale may also need to be changed when the probe is on a narrow stenosis, causing the peak velocity to increase by four times. As soon as a stable measurement is reached, or whenever the surgeon says so, save the measurement by pressing the "Save" button.


5.5.4 Transit time flow measurements

The VeriQ™ displays the flow curve at a defaulted 5-seconds sweep rate upon probe connection. Look for a green Acoustical Contact Indicator, which indicates appropriate contact between the probe and the vessel. Inform the surgeon of the mean flow rate, the pulsatility index (PI) and the diastolic filling (DF). As soon as a stable measurement is reached, or whenever the surgeon says so, save the measurement by pressing the "Save" button.

5.5.5 Optimizing the Doppler and flow display

Change of vertical scale

If the Doppler spectrum or transit time flow curve appear very small, or exceed the height of the display, the vertical scale should be changed.




1. Press on the top left part of the scale to increase the scale. This will decrease the height of the displayed signal.
2. Press on the bottom left part of the scale to decrease the scale. This will increase the height of the displayed signal.

Adjusting baseline position

The baseline position can be changed by dragging the baseline up or down. The Doppler baseline can only be adjusted when recording "live", while the flow and AUX channels can be changed in "Pause", "Edit" and live recording modes.

Change of time sweep

Press either one of the icons on the top right side to change the horizontal time scale.



1. Press (1) to increase the time sweep.
2. Press (2) to decrease the time sweep.

5.5.6 Saving Measurements


When saving, the last one minute of the recorded trace will be stored. During a search for vessels, only the velocity trace of the located vessel is necessary to save. When detecting a stenosis, it will be important to save both the reference velocity trace and the trace from the position of the stenosis. A transit time flow measurement can be saved as soon as a stable measurement is reached.



Press the "Save" button on the icon menu to store recorded traces.

5.5.7 Memory scroll and trace review

The measurements need to be paused before memory scroll can be performed.



1. Press the "Pause" icon on the icon menu (1) to halt the measurements.
2. Drag a finger across the display to review recorded data.
 - a. Drag a finger from left to right to move back in time.
 - b. Drag the finger from right to left to move forward in time.

Allow some time for the system to recognize the finger position before scrolling. Press the "Save" icon (2) to store the current position

Memory scroll can be used in the "Archive" screen or "Paused" mode when editing curves.

5.5.8 Review of recorded data


All stored measurement data can be reviewed. A thumbnail of each stored data-set is presented on the screen, and the number of thumbnails visible on a single page can be adjusted.

System settings description

1. Press the "Data screen" tab (1) to enter the explorer screen.
2. Use the scroll bar (2) to view any thumbnails may not be visible on the screen.
3. Select the wanted thumbnail for additional information

by pressing on the thumbnail (a check will appear on the selected thumbnail).

4. Press the "Edit" button (3) to review and/or edit the selected thumbnail.
5. Press the "Printer" button (5) to print the patient report with the selected thumbnails.
6. Press the "Export" button (6) to export the report to a desired external memory device. The report will be exported as a PDF file, with the same data as the printout.
7. Press the "Delete" (7) button to delete the selected thumbnail(s).



Note

After editing a curve, changes can be made permanent by pressing the "Save" icon. When multiple curves are selected for editing, pressing the arrows on the bottom of the screen will bring the next selected file into edit mode.

Note

All thumbnails can be selected by pressing the box containing the date and time of operation.



Press the "Back" icon to re-enter the "Archive" screen when in editing mode. For the more-advanced user, an "Advanced" system mode can be chosen.

5.6 Tasks during operation by surgeon

5.6.1 Doppler velocity measurements


DP ✓ 45°

To get an accurate velocity measurement, the Doppler probe should be held perpendicular to the vessel being measured, as its transmitting crystal is fixed at a 45° from the probe's surface. Movements of the probe, such as those caused by a beating heart, will cause a strong low-frequency noise on the Doppler spectrum. This can be avoided by applying a mechanical stabilizer, or for smaller movements, by increasing the low velocity filter in the system.

The Doppler probe should be held with the angled crystal pointing against the flow direction, as shown on the illustration above. If the probe is held the opposite way, the velocity curve will appear below the baseline.

The operator can change the filter settings by pressing the Doppler probe identity button.

For correct data management, it is important to inform the system operator about the measured vessel name and the probe position on the vessel (Ex. LAD Proximal, LAD Stenosis or LAD Distal) and to inform the operator when to save the recorded Doppler spectrum.




5.6.2 Flow measurements

As with the Doppler measurements, the operator needs to be informed which type of conduit, connection, and target vessel is being measured. When the flow probe is placed onto a conduit, the system will generate a simulated Doppler sound, indicating that flow has been detected.

The volume of this sound is adjustable, using on-screen buttons. The operator should confirm that ACI (Acoustical Coupling Index) value is > 50%, with a green background, ensuring that measurements will be valid. If the ACI is not green, consider using sterile ultrasound couplant, spraying saline on the probe head, skeletonizing the vessel or using a different probe size. If there is a gap between the vessel wall and the probe surface, the probe is too large and a smaller probe should be used. If parts of the vessel are squeezed out of the probe, the probe is too small and a larger probe should be used.

To ensure accurate flow measurement, it is important that the probe is placed perpendicular to a straight segment of the measured vessel. All parts of the vessel should be inside the probe.



Recommended probe placement on straight vessel segments.

5.7 Post-operative procedures

5.7.1 Retrieve patient from database

1. Go to the Home screen by pressing the "Home" icon
2. Press the "Search" button
3. Specify search criteria
4. Press "Perform search"
5. Select a patient by clicking on his/her name and then "Select Patient" at the bottom of the page. To export the patient data, click on the patient's corresponding blank box under the heading "Select".
6. Click on "Export"
7. Press the "Data" icon to view the stored data for the selected patient

5.7.2 Cleaning and disinfection of probes

It is important to clean the probe immediately after the operation, to prevent blood from drying on the probe. During cleaning, the probe must be handled carefully, avoiding stretching of the cable or bending of the neck or steel reflector on the probe. The probe can be cleaned in lukewarm water using a mild soap-solution.

If the probe needs to be disinfected before handling, it is important that only solutions tested and validated by Medistim ASA are used on the probe. An updated list of tested solutions may be found on our web site, www.medistim.com

5.7.3 Sterilization

The different Medistim probe series must be sterilized using specific methods for each series. See Appendix B.2 for an overview of the probe cleaning manuals available for each probe series.


Note

The probes have been limited to 50 sterilization cycles / usages. When a probe reaches 50 usages it will be deactivated.

5.7.4 Reporting

Following the operation, the surgeon should review the stored data in the "Archive" screen, and select the measurements that best represent the performed surgical procedure.

The selected measurements will be marked with a green check-mark, indicating that they are selected for reporting. The report can either be printed directly by pressing the "Print report button" or stored as a PDF file locally or on a removable USB stick.



The report layout can be changed in the "Report layout" pop-up.

Report Settings

Report Layout



☒ Show report layout selector when generating a report.













Ok

Measurements that have been reported will be flagged with a green check mark and the text "Qualified". These "Qualified" measurements should not be changed as they are now regarded as valid documentation for the performed procedure.




5.8 Icon description


The following overview gives a description of all icons used in the different system screens.


	Home data screen - Screen to record patient data, edit patient information and retrieve saved information		Flow wizard - Displays standard flow curves when commonly-used vessel names are chosen.
	Measure screen - Screen for Doppler and flow measurements, identifying vessel names and saving measurements		Invert - If the flow curves goes below the baseline, pressing this button will invert the curve.
	Imaging screen - Screen for Imaging measurements. Do 2D, CFM or PW-Doppler, identify vessels and save measurements. Only available on VeriQ C™.		Export report - Selected measurements are exported to an external memory device in PDF format.
	Data screen - Edit, delete, print and export saved measurements.		Print - All selected measurements will be printed.
	Help screen - Video tutorials and the electronic User Manual for the VeriQ™ System can be found here.		Delete - Will delete selected item(s).
	System Settings - Customize system setup here.		Select / Edit - View or edit selected measurement curves.
	Time Compression - Decrease time sweep for all traces in steps. Down to 1 second.		Display 4 measurements at a time. Displays measurements with large thumbnails in a 2 by 2 arrangement.
	Time expansion - Increase time sweep for all traces in steps. Up to 60 seconds.		Display 9 measurements at a time. Displays measurements with large thumbnails in a 3 by 3 arrangement.
	Save Data - Save recorded data to the patient database.		Display 16 measurements at a time. Displays measurements with small thumbnails in a 4 by 4 arrangement.
	Measurement pause - Press to pause measurement. The icon will change to the play symbol.		Thumbnails - Go to explorer screen. Press this icon to return to the explorer screen when in editing mode.
	Measurement play - Press to resume measurement. The icon will change to the pause symbol.		Screen shot or video export - The selected screen shot or a video of the ultrasound recording can be exported to an external memory device.
	Derived curves - Choose combination of flow and pressure channels for derived curves.		Move up one level when browsing the file system.
	Trend setup - choose number of channels and set sweep time for trend recordings.		Create new folder.
	Trend / live - In trend measurement mode, this button will allow the user to switch between trend and live mode measurements.		Video tutorials - Press this icon to access video tutorials on the VeriQ system.
	Measure and Analyze - Displays calculated values and measurement tools.		User manual - Press this icon to access the electronic user manual.
	Volume control - The lowest position will turn sound off.		Increase Doppler gain - The spectrum will become darker with the possibility of increased noise. Used with the Medistim Doppler probe.
	Sound on / off		


Q1 4mm ACI 100 %


Probe properties - Probe information, filter settings, liquid type, temperature and sound, velocity settings, Measure and Analyze, index selection and ACI.


 **Decrease Doppler gain** - The spectrum will become more transparent with less noise. Used with the Medistim Doppler probe.


 **Increase sample volume** - Increases the sample volume of the Doppler by 1mm.


 **Decrease sample volume** - Decreases the sample volume of the Doppler by 1mm.

 **Doppler reference speed** - This button will mark the current maximum speed and use it for stenosis grade calculations.

 **Search Archive** - Triggers a search operation and searches for patient(s) in the internal archive.


 **New patient** - Triggers the creation of a new patient in the internal archive and will open a “new patient” screen.

 **Edit patient info** - Edit or add information to a patient file.

 **Comments** - Opens a comments field for the current measurement set.

5.9 System settings description

Under the tools icon there is a comprehensive system settings menu, controlling system configurations and feature not frequently used.



The system settings menu contains the following settings:

1. Name and address of hospital can be entered/ edited and will be recorded on the patient report
2. Surgeon setup, with drop-down menu to enter/ select name of surgeon
3. The vessel name setup allows editing or adapting factory pre-stored vessel names
4. Report settings allows selecting desired number of curves per page in the patient report

5. Probe report indicates the usage status of each probe being used on the system
6. The manage projects setting allows creating projects names for sub-groups of patients
7. The advanced system functionality setting allows for management of patient database and other rarely used system operations

There are 4 pull down menus:

1. **System profile** for selecting the AUX channel used for the ECG signal
2. **ECG channel** selecting the AUX channel used for the ECG signal
3. **Measurement unit** allows the displayed unit to be changed from metric to US
4. **System mode** allows hiding the most advanced functions for easier operation

There are also 5 optional selections:

1. Show scale buttons displays the buttons to increase and decrease vertical scale
2. Use automatic vessel pop up will prompt use for vessel name each time a probe is connected
3. Enable Auto Invert automatically invert the flow curve when mean flow is below baseline
4. Show Button Help Text explains the function of each screen icon.

5.10 Other adjustments

For adjustments such as setting the time and date, name of surgeons, editing the vessel names, deselecting automatic vessel pop-up and adjusting the number of thumbnails visible in the “Archive” screen.

There are two system modes: advanced and clinical. The advanced mode contains options such as derived curves and trend recordings. The clinical mode is tailored for common cardiac and vascular surgeries, and only the most frequently used options are available. The advanced and clinical user modes are activated in the system settings menu.

5.11.1 Advanced Doppler mode

If the probe controls are not visible in the Doppler display there are two options:

1. Press the button reading “DP V 45°” and adjust the settings in the screen that appears.
2. Press the button reading “DP V 45°”, select “Show Probe Controls in Curve View” and adjust the gain in the Doppler display. The selection will be active until it is de-selected.

Optimizing display gain

When the Doppler spectrum appears faint or saturated, the gain can be adjusted.




1. Press the "+" button (1) in the lower right part of the Doppler display to increase the gain if the spectrum appears faint.
2. Press the "-" button (2) in the lower right part of Doppler display to decrease the gain if the spectrum appears noisy.

Change of sample volume

By default the system measures all velocities at depths between 2 mm and 8 mm. The field of measurement and the sample volume can be changed.

To change sample volume

1. Press the "+" button (3) in the Doppler display to increase the sample volume
2. Press the "-" button (4) in the Doppler display to decrease the sample volume.
3. To increase or decrease depth, press the top or bottom of the scale respectively.



As shown in the figure above, the number on the right indicates the depth of the measurement in mm and the number on the left indicates the sample volume in mm.

Note


The depth is measured from the probe surface, and indicates the real depth to the center of the sample volume. However, since the beam is transmitted at a 45° angle, the actual measured sample is situated in front or behind the probe depending on probe direction.

Change the filter settings

Probe movements may introduce noise close to the baseline, especially when measuring on a beating heart. Increasing the filter setting can reduce the amount of noise. Press the "DP V 45°" button and select a higher filter value under "HP Filter Frequency". The factory default is 100hz.

5.10.1 Tuning Doppler noise

The Doppler measurement is sensitive to noise from the instruments nearby, possibly causing the calculated values from the measurement to be inaccurate. Therefore this calibration should be done in an environment resembling the working environment of the VeriQ™ in question. For the purpose of tuning of the noise reduction, the probe can be placed in a cup of water.




Turn on the VeriQ™ and insert a Doppler probe. Enter the Doppler settings panel.

Press the "Advanced" button. The "Noise Suppression setting" should now be visible.

The slider indicates the level of noise reduction to be applied to the measurement. A negative value will make the measurement more sensitive to weak signals and noise, and a large positive value will make it less sensitive. Setting the noise reduction too high may cause problems calculating measurement parameters, as the measured signal is interpreted mostly as noise. Similarly, setting the value too low will result in much noise being incorporated into the calculations. Adjust the slider so that the probe stops measuring the environmental noise.

The noise reduction value can at any time be readjusted if the value initially chosen is causing problems or the environment changes.



5.11 Data Management

All measurements are stored in a database associated with the patient's record. The capacity of the database is typically about 40,000 records, depending on the number of measurements stored per patient.

5.11.1 Search function

With the option of multiple search criteria, any desired record is easily found. To search for a particular patient, operation date or surgeon, go to the "Home" data screen and select "Search". The screen "Search Criteria" appears. Search for a particular patient by entering parts of the patient name, the patient ID number, the birth date, the gender, the operation date or the surgeon. Press "Perform Search". Choose the desired patient record from the list that appears by pressing on the name - a green check appears - and pressing "Select Patient".

Note

The project can be used to identify patients part of a study making search easier.

5.11.2 Edit, print, delete and export measurements

To access a patient's measurements, make sure the name of the patient is selected, and go to the "Archive" screen. View the desired measurements by pressing on the curves.

Press the "Edit" button and edit the curves by either scrolling the flow curve backwards, dragging the baseline up or down, changing the scale of the measurements, or editing the vessel name.

The scale and the position of the baseline cannot be changed on Doppler measurements. Move between the selected measurements by pressing the arrows visible at the bottom of the screen. To store the edited curves, press the "Save" button. The old measurement will then be overwritten.

To print measurements, select the measurements to print by pressing on the curve and pressing the "Print" button. To change the number of measurements visible on a printed page, change the setting under "Report Settings" in the system settings.

To delete, select the measurements by pressing on the curve and then press the "Delete" button. The question "Delete selection?" will appear, press "OK".

To export measurements in a PDF format, select the desired measurements and press the "Export report" button. The screen "Specify File-name and Folder" will appear. Choose the desired folder into which the file is to be exported, specify the name of the file and press "OK". The report can be exported to the hard disk, an external USB or LAN storage device.

The full screen can be exported as a picture file for inserts into presentations or other documents. Press the screen capture button for the edit screen and follow the dialog box to store the picture file on a removable media.

5.11.3 Export and import data

It is possible to export either the entire patient database, or portions of it, for backup purposes. This can be done by using an external USB or LAN storage device. To export single patients, choose Advance Mode under System settings to gain access to the export function. Go to the "Patient data" screen and obtain the full list of patients by pressing "Search" and then "Perform Search".

Under the heading "Export", select the desired patient data and press "Export". Multiple patients can be selected. Choose the desired folder under "Specify Filename and Folder" and press "OK".

On the advanced systems settings page there is an option to export various types of data, including the entire patient database. The data available for export is:

- Entire patient database
- Language files
- System settings
- GUI configuration
- Probe usage information

To import data, enter the "Advanced functionality" menu found in the system settings and select "Import". Select the desired files and press "OK".

Note



The raw patient data files can only be read by VeriQ™ or Medis-tim workstation program.

5.11.4 Safe removal of storage devices

Removable storage devices must be prepared for safe removal before disconnecting from the system.

If a removable storage device is disconnected while the system is transferring or saving data there is a risk of data loss or data corruption.


The option for safe removal will appear automatically but can also be selected under “System Settings”. Once a removable storage device has been prepared for safe removal, it must be disconnected and then reconnected in order to be detected by the system for further use.



5.12 Derived traces

In the advanced mode it is possible to calculate derived traces. This is done by adding, subtracting or dividing the traces of two inputs. Resistance can be calculated by dividing pressure by flow. When measuring the carotid bifurcation, the flow in the Carotis Externa and the Carotis Interna can be derived into one total flow amount.

For enabling derived traces, enter the “Measure screen”, insert the desired flow and pressure probes, and press the button “Derived curves”. Choose the desired calculations of traces and press “OK”.




5.12.1 Fast Fourier Transformation (FFT)

Included in the M&A software is a capability of doing Fourier analysis (frequency analysis). The analysis can be performed on all flow, pressure and AUX curves, including derived curves. By using this capability it is possible to analyze the frequency content of the physiological signals, changing the x-axis into a frequency axis with the unit Hz. The signal is converted into different columns along the x-axis, showing where the different frequency components are located. The height of the different columns is proportional to the intensity of the different frequency components. The highest frequency component is normalized to 1 on the display.

It is recommended that a long time signal, 30 seconds or longer, is used when performing a Fourier analysis of a physiological signal. The Fourier analysis is done just by pressing the Derived Curves Button.

In the pop-up menu, select the channel to be FFT analyzed. The result of a Fourier analysis will be shown in a new window at the bottom of the screen.

The Fourier analysis can be stored and printed. M&A calculations can also be performed on an analyzed signal. By moving the cursors C1 and C2 along the x-axis, frequency and intensity of the different frequency component will be updated in the M&A window.



5.13 Trend measurements

Recording the flow trend over longer periods of time is another option available in the advanced mode. For example when continuously measuring aortic output, trend measurement will display the flow development over time.

For enabling trend measurements, enter the "Measure screen", insert the desired flow probe, and press the button "Trend setup". Choose the correct flow channel under "Available channels" and the desired sweep length under "Display sweep length". Press "OK".


In trend mode, switch between the trend and the live mode by pressing the "Switch between trend and live mode" button. To stop the trend logging, press the "Trend setup" button and choose "Stop logging". Answer "OK" to the question "Stop trend data collection?" The trend measurement will be automatically saved.

5.13.1 Save, review and edit trend curves

In the trend mode the "Save" function works differently from the live mode. Pressing "Save" will save the trend measurement and the event at that exact time. The purpose of this is that both the trend as well as the timing of events, e.g. injection of medication, can be reviewed.

Pressing "Save" will make the screen "Event description" appear. Fill in a description of the event and press "Done". Repeat this until all the events are saved and the desired sweep time. The trend measurement is automatically saved every minute, and is stored as a whole after logging.

To review or edit the trend curves, go to the "Archive" screen. Select the measurements to review or edit, and press the "Edit" button.




Note

For the best possible result on the FFT it is important to have stable measurements for 30 seconds or more.

Events are highlighted in pink on the trend measurement curve. To view description of events, press the "Events" button at the bottom of the screen. This text may be edited by clicking in the text box. To return to the trend measurement, press "Trend". The events are also stored as separate flow curves, which may be reviewed, edited and printed like any flow measurement. To print trend measurements, select the desired measurements and press "Print".

5.14 Measurement and Analysis (M&A)


Pressing the button labeled "Calculate" with an icon symbolizing a calculator enters the Measurement and analysis mode. In this mode various different calculations and indexes related to the measurement curves are shown on a pane to the left of the curves.



Pressing the probe button brings up the property sheet for the measurement. Here it is possible to select or deselect the different measurements.

Note

Not all indexes/measurements are visible in live mode.



In Edit mode all the indexes are visible and by moving the cursors C1 and C2 the corresponding values are updated. The values shown for all the indexes, except for C1 and C2, are based on the interval between the two cursors. For instance, the "Max" index displays the highest value measured in the interval between the cursors.

6. TROUBLESHOOTING

6.1 Service of the VeriQ™ System

For service instructions, consult the VeriQ™ System service manual or contact the Medistim service department directly.

6.2 Basic troubleshooting

The following are some basic tips in the case that the VeriQ™ System is not operating properly. Contact the hospital technician if these tips do not resolve the issue.

- Ensure that the system is connected to a functioning power outlet and that the system is switched on.
- Verify that all probe connections are properly inserted at the junction point.
- Visually inspect probes for damage and replace with a new probe if necessary.
- Confirm that all settings are done in accordance with the instructions in this User Manual.
- Verify that connections are made in accordance with the wiring diagram.

Note

Refer to the system's part number and serial number when reporting any issues.

6.2.1 Troubleshooting users guide

Symptom	Possible failure mode (investigate in listed order)	Resolution
MONITOR		
No picture on monitor (power button green)	1. Monitor powered off	Use power switch on the monitor if needed, check OSD settings, (refer to Monitor's user's Guide)
	2. DVI-D Cable or Monitor power cable not properly attached	Connect the cables to the monitor
SYSTEM OPERATION		
System does not start when pressing power button	System not connected to mains	Check that power plug is properly connected to a powered grounded outlet
CONNECTING PROBES		
Probes / Aux or Pressure contacts cannot be inserted into port	Pinhole(s) of port obstructed	Use other port if available. Report to technician. Until resolved, block failing port with tape.
Flow/Doppler Probe not detected when inserted or Pressure/Aux signals not detected or System message "probe failed" or 0% ACI.	1. Probe connector not fully inserted into port	Check that probe connector is fully inserted
	2. Broken/bent pins inside the probe connector	Visually inspect the connector
	3. Other probe damage	Try using probe on other port of same type. If the problem persists, dispose and replace with a new probe.
PRINTING		
Poor printing quality	Out of Ink	Replace ink cartridge. (Refer to printers User's guide)
No print is generated	1. No paper loaded	Load paper. (Refer to printer User's guide)
	2. Printer has accidentally been turned off	Switch the printer on (Refer to printer User's guide)
	3. Printer power or USB cables unplugged	Check that the cables are correctly attached. (Refer to printer User's guide)
TOUCH MONITOR RESPONSE		
Incorrect function evoked when pressing the screen	Touch screen dirty	Clean the screen
No response when pressing the screen	Touch USB-cable disconnected from monitor	Reconnect cable and reboot the system.

7. MAINTENANCE

The purpose of this section is to give the user applicable information regarding maintenance and warranty requirements. Additionally, user-replaceable parts and miscellaneous repair procedures are presented.

7.1 Probe statistics

The various probes which have been connected to the system are automatically identified, and serial-number specific probe data is stored. Probes are limited to 50 uses and can be removed and re-connected as many times as desired during a five-hour period, with only one usage being counted. If the system is turned off and restarted again during these five hours, the probes will count a new usage when re-inserted.

The system automatically generates a message when there are 5 uses left. On the last use, the system gives a probe expiration warning. When connecting a probe which has been used 50 times, the system will reject the probe.

The system keeps track of the usage of all probes. To obtain probe usage statistics, use the following procedure:

1. Go to System settings. Push the "Probe Report" button.
2. Information for all probes appears on the screen. Different filters can be applied

The probes are all identified by a discrete serial number. The probe list can be sorted by pushing the header row buttons. The toggle buttons in the "Select" column can be used to select probes for printing or deletion. Pushing the "Select" button in the header row will select/deselect all probes in the list.

Pushing the "Delete selected" button will delete all currently selected probes from the list. (The deleted probes may still be listed by using the "Deleted probes" filter).

The "Print selected" button prints a probe report containing the selected probes. The following filters are available:

Patient 01/23/14 09:50:59 23. Jan 10:01 MEDISTIM

Probe Report

Filter: All probes

Serial number ↑	Probe size	Probe type	Last inserted	Uses	Remaining uses	Select
EL10001500057	STI-15MHz		23.01.2014	148	Unrestricted	<input type="radio"/>
EL10001500023	STI-15MHz		23.01.2014	512	Unrestricted	<input type="radio"/>
PQ000220131	2,0 mm QuickFit		23.01.2014	13	37	<input type="radio"/>
PQ00032666	3,0 mm QuickFit		23.01.2014	50	0	<input checked="" type="radio"/>
PB0003240434	3,0 mm PeriVascular		23.01.2014	13	37	<input type="radio"/>
PQ0003208216	3,0 mm QuickFit		23.01.2014	9	41	<input type="radio"/>

Delete selected Report selected OK

Button

"All Probes"

"Safety check required"

"Deleted probes"

Explanation

Lists all (non-deleted) probes used with the system .

Lists all imaging probes that require regular safety checks (used more than 100 times).

Lists all probes that have been deleted from the list.

7.2 Periodic maintenance

- The system surface can be cleaned with hospital grade disinfectants. Avoid excessive use of water or liquids on the system.
- Verify that the main power cable and connector are not damaged.
- Check that the cable attachments to the monitor and the electronic unit are properly mounted.

7.3 Storage

The probes should be kept in their container when not in use. The system and probes can be stored under the following conditions:

- Temperature between -25 °C to 70 °C (-13 °F to 158 °F)
- Humidity between 20 % - 95 % RH, non condensing

When the system or probes are moved between locations with different temperature and humidity, please allow sufficient time for acclimatization to preventing condensation.

7.4 Packing and transport


For safe transportation, please use only original packing materials. Note that a special shipping container, specifically designed to protect the electronic unit and monitor assembly, is available from your local distributor.

7.5 Disposal instructions

The system contains various plastic materials, electronic components and metals. Please follow local rules and regulation when disposing. Please follow the hospitals safety regulations for equipment that has been in contact with blood.

8. TECHNICAL SPECIFICATIONS

8.1 VeriQ™ System overview



SOCKET FOR SYSTEM
PRINTER ONLY

POWER INLET WITH
FUSE HOLDER

SOCKET FOR FRONT
END UNIT ONLY (VERIQ
XXXXC SERIES)

SOCKET FOR
SYSTEM UNIT (PC)
ONLY



DO NOT CONNECT ANY
OTHER TYPE OF
EQUIPMENT TO THIS
SOCKET.

8.2 Technical specifications

PHYSICAL PROPERTIES		
Dimensions System complete		58 x 140 x 58m (WxHxD)
Monitor:		38 x 32 x 6 cm (WxHxD)
System weight		VeriQ™: 56.2 kg VeriQ C™: 58.6 kg
Ingress protection		
System:		IPX0
Transit-time probes:		IPX7
PV probe:		IPX7
Imaging probe EL100015:		IPX7
Imaging probe EL100010:		IPX1
ELECTRICAL SPECIFICATIONS		
Supply voltage rating		220 V- 240 V ~ , 50 - 60 Hz, single-phase 100 V - 120 V ~ , 50 - 60 Hz, single-phase Max 500VA
Supply power rating		230 V models: 2 x T5AH
Fuses		115V models: 2 x T5AH
Patient safety		
Patient applied parts		Type CF, defibrillator proof
Medical safety		Compliant with medical standard IEC 60601-1 Class I, protectively earthed.
CLEANING, DISINFECTION AND STERILIZATION		
Cleaning		The system surface can be cleaned with a soft damp cloth or if necessary a mild non-alcoholic soap. Avoid excessive use of water or liquids on the system. The touch screen area can be cleaned with anti static cleaning tissues. For cleaning of probes before sterilization consult the “Handling, Cleaning, Disinfection, and Sterilization of the Medistim Probes” packaging leaflet for each probe type.
Disinfection		Approved disinfectants are specified in the “Handling, Cleaning, Disinfection, and Sterilization of the Medistim Probes” packaging leaflet for each probe type.
<i>Note that any approved disinfectants are recommended because of their chemical compatibility with product materials, not their biological effectiveness. For the biological effectiveness of a disinfectant, refer to the guidelines and recommendations of the disinfectant manufacturer.</i>		
Sterilization		The probes can be sterilized using low temperature gas-plasma. Some probe types also supports steam autoclave. Consult the “Handling, Cleaning, Disinfection, and Sterilization of the Medistim Probes” packaging leaflet for each probe type.
ULTRASONIC PARAMETERS		
Reporting and compliance model	FDA:	IEC 60601-2-37:2007 Track 1 (VeriQ, and VeriQ C™ TTFM and doppler probes) Track 3 (VeriQ C™, imaging probe)
Maximum acoustic levels		
Transit-Time (VeriQ™, VeriQ C™)		MI < 0.02 Ispta.3 < 0.1mW/cm2
PW-Doppler (VeriQ)		MI < 0.05 Ispta.3 < 20 mW/cm2
B-Mode (2D) (VeriQ C™)		MI < 1.1, TIS < 0.9
CFM (VeriQ C™)		MI < 1.8, TIS < 5.2
PW-Doppler (VeriQ C™)		MI < 1.9, TIS < 4.3
Acoustic exposure parameters		MI, TIS (VeriQ C™ only)

Ultrasonic parameters

Transit Time	Transducer frequency: Excitation: Repetition rate (PRF)	500 kHz - 3.5 MHz Burst of 16 - 25 waves 0.6 kHz - 2 kHz	
PW-Doppler (VeriQ)	Transducer frequency: Repetition Rate (PRF) :	7.5 MHz 4 kHz - 32 kHz	
Imaging (VeriQ C™)	No. of physical channels: Transducer frequency: Ultrasound bandwidth: Maximum frame-rate: Supported modalities:	128 11 – 18 MHz 20 MHz 50 fps 2D (B-Mode)	CFM-Velocity Doppler CFM-Power Doppler PW-Doppler
Supported probe geometries:		Linear, phased array	

REGULATORY COMPLIANCE

Sales Classification	Rx-Only. Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician		
Device Category	Medical ultrasonic imaging and volume flow-meter system with probes		
Quality System	Medistim ASA is ISO13485:2003 certified		
Provisions Imaging and volume flow meter probes	Risk Class III in Annex II of the EC-Directive 93/42/EEC (MDD) flow-meter probes: as amended by Directive 2007/47/EEC (CE-mark)		
Provisions VeriQ™ and VeriQ C™ Systems	Risk Class IIb in Annex II of the EC-Directive 93/42/EEC (MDD) as amended by Directive 2007/47/EEC with the exemption of section 4 (CE-mark).		

SYSTEM COMPONENTS

System Monitor	Type: Touch screen: Native resolution: Power:	19" TFT LCD, IEC60601-1 compliant Surface-wave type, requires fingertip operation only, sharp pointing devices will not work. 1280 x 1024 Pixels 12V DC, 50W
CPU	Intel Core2 dual core	
Memory	4GB	
System Storage	System hard disk >= 500GB (one additional HDD for data storage is optional)	
Smart card reader	ISO/IEC 7816 compatible	
Ports	2 x USB 2.0	
Sound	Stereo speakers for directional Doppler sound	
Printer	InkJet printer, colour, USB 2.0 compatible	

MATERIALS

Component	Material	Coating/Surface
Base frame	Molded aluminium	Powder paint
Columns	Extruded aluminium	Powder paint
Power inlet module	Sheet metal	Powder paint
Front cover	Sheet metal	Powder paint
Top cover	Molded PUR	Paint
Back cover	Thermoformed ABS	Structured (VeriQ™) Paint (VeriQ C™)
Other structural brackets	Sheet metal	Aluzinc
Handles	Thermoplastic	

STANDARDS

Safety	IEC 60601-1 (ed.3)
Usability	IEC 60601-1-6 (ed.3) / IEC 62366 (ed.1)
Ultrasound	IEC 60601-2-37 (ed.2)
EMC	EN 60601-1-2 (ed.3)
Biocompatibility	ISO10993-1:2009
RoHS	EN50581 (ed.1)

APPROVALS



The Medistim system and probes are approved according to UL 60601-1 & CSA C22.2 No. 601.1 with MET listing number E113031.



The Medistim VeriQ™ Systems and probes are CE approved.

SYSTEM SOFTWARE

Operating system	WindowsXP Professional SP 3
Application software	VeriQ™ System Software v3.31 or later (VeriQ™, VeriQ C™)

NETWORK AND CONNECTIVITY

Network connection	1000BASE-T (IEEE 802.3ab)
Protocol(s)	DICOM (optional, refer to DICOM Conformance Statement)

FLOW CHANNELS

Frequency response	DC to 100 Hz
Acoustical Coupling Indicator (ACI)	Probes are calibrated to 100% in water Green >50 % Yellow 30 % - 50 % Orange 10 % - 30 % Red <10 % (signal are forced to zero)Screen output
Pulsatile and mean flow curve	Mean flow value, pulsatility index and diastolic filling index
Measure and Analysis	Maximum, minimum value, insufficiency rate, flow integration, maximum derivative and cursor readings
Sample rate	0.3– 1kHz depending on probe
User settings	Noise filter, mean averaging time, scale, sweep-length and signal inversion
2D (B-Mode) Image (VeriQ C™)	
Screen output	Image, highest Tx frequency, frame rate, exposure indices
User settings	Gain, persistence, frame-rate/resolution priority, depth of view, focal point selection.
Color Doppler (VeriQ C™)	
Modes	Velocity and power Doppler
Pulse Repetition Frequency	0.5-6 kHz
Screen output	Colour map, frame-rate, exposure indices, 2D and Doppler Tx frequency
User settings	Colour gain, Region Of Interest (ROI), Beam angle, PRF
Pulsed Wave (PW)- Doppler	
Wall filter	50-1500 Hz in steps
Pulse Repetition Frequency	VeriQ C™: 3-16KHz VeriQ™: 4-32kHz
Screen output	Velocity spectrum, maximum spectral velocity, mean spectral velocity, average max spectral velocity index percentage velocity difference
Measure and analysis	Maximum and minimum value for maximum spectral velocity and mean spectral velocity, velocity integration, maximum derivative and cursor readings.
User settings	Wall filter, PRF, sample volume, gate positioning and receive gain, curve display selection

PRESSURE CHANNELS

Transducer sensitivity	5 μ V/V/mmHg
Range	360 mmHg
Accuracy	$\pm 2\%$
Resolution of indexes	1 mmHg
Frequency response	DC to 100 Hz
Sample rate	250 Hz
Screen output	Pulsatile and mean pressure curve Systolic, diastolic and mean pressure value
Measure and Analysis	Maximum, minimum value, maximum derivative and cursor readings
User settings	Zero-point, gain and noise filter

AUXILIARY CHANNELS

Signal input ranges	± 10 mV ± 100 mV ± 500 mV ± 4 V
Sample rate	250 Hz
User configuration	Input range, Offset, scale, axis, unit, noise filter
Settings storage and handling	Fully profile based. Systems support 10 named profiles per channel and selection of default profile at power up
Screen output	Pulsatile and mean curve, mean value
Measure and Analysis	Maximum, minimum value, maximum derivative and cursor readings

GENERAL SETTINGS FOR ALL CURVES

Noise filter	5, 10, 20, 30, 50, or 100 Hz
Measure & Analysis Indices	User selectable indices to be displayed
Mean averaging time	0.1 - 10 seconds averaging time, selectable Factory setting is 7 seconds
Sweep length	1-60 seconds
Max measurement save length	60 seconds for pulsatile curves. Trend measurements are limited only by procedure length.

ENVIRONMENTAL CONDITIONS

Storage / transport	Temperature Humidity Atmospheric Pressure	-25°C to 70°C (-13°F to 158 °F) 20% RH to 95% RH (non-condensing) 101 kPa +10% to -20%
Operation	Temperature Humidity Atmospheric Pressure	0°C to 40°C (32°F to 104°F) 20% RH to 80% RH 101 kPa +10% to -20%
Transport	Use original packaging	

APPENDIX A. PROBE APPLICATION OVERVIEW

A.1 Probe application and size guide

Cardiac Surgery

Surgery/Vessel	Vessel Size	Recommended Medistim probes
Saphenous Vein Graft	3-5 mm	3mm, 4mm and 5mm Medistim TTFM probes (PS and PQ probe series)
Internal Mammary Artery	2-3 mm	2mm and 3mm Medistim TTFM probes (PS and PQ probe series)
Radial Artery	2-4 mm	2mm, 3mm and 4mm Medistim TTFM probes (PS and PQ probe series)
Ascending Aorta	20-35 mm	Please contact your sales representative or Medistim directly for information
Pulmonary Artery	20-35 mm	Please contact your sales representative or Medistim directly for information
Pediatric Surgery	6-20 mm	Please contact your sales representative or Medistim directly for information

All probes are available with or without handle. See section A.2 for a full list of available probes.



Transplant Surgery

Surgery/Vessel	Vessel Size	Recommended Medistim probes
Hepatic Artery	5-8 mm	5mm, 6mm, 7mm and 8mm Medistim TTFM probes (PV, PS, or PQ probe series)
Portal Vein	10-14 mm	10mm, 12mm, and 14mm Medistim TTFM probes (PV probe series)
Renal Artery	4-6 mm	4mm, 5mm and 6mm Medistim TTFM probes (PV, PS, or PQ probe series)
Renal Vein	8-11 mm	8mm, 10mm, and 12mm Medistim TTFM probes (PV probe series)
Common Iliac Artery	6-8 mm	6mm, 7mm, and 8mm Medistim TTFM probes (PV or PS probe series)

All probes are available with or without handle. See section A.2 for a full list of available probes.

Vascular Surgery

Surgery/Vessel	Vessel Size	Recommended Medistim probes
Common Carotid Artery	6-9 mm	6mm, 7mm, 8mm and 10 mm Medistim TTFM probes (PV and PS probe series)
Internal Carotid Artery	4-6 mm	4mm, 5mm and 6mm Medistim TTFM probes (PS and PQ probe series)
External Carotid Artery	4-6 mm	4mm, 5mm and 6mm Medistim TTFM probes (PS and PQ probe series)
Common Femoral Artery	8-11 mm	8mm, 10mm and 12mm Medistim TTFM probes (PV probe series)
Popliteal Artery	3-6 mm	3mm, 4mm, 5mm and 6mm Medistim TTFM probes (PS and PQ probe series)
Tibial Artery	3-4 mm	3mm and 4mm Medistim TTFM probes (PS and PQ probe series)
Saphenous Vein	4-6 mm	4mm and 5mm Medistim TTFM probes (PS or PQ probe series)
Renal Bypass	4-6 mm	4mm, 5mm and 6mm Medistim TTFM probes (PS or PQ probe series)
Descending Artery	20-30 mm	Please contact your sales representative or Medistim directly for information
Radial Artery	2-3 mm	2mm and 3mm Medistim TTFM probes (PS or PQ probe series)
Brachial Artery	3-4 mm	3mm and 4mm Medistim TTFM probes (PS or PQ probe series)



A.2 List of available Medistim probes

Medistim QuickFit™ TTFM probes - PS probe series

Probe name	Probe sizes	Part numbers*
QuickFit™ TTFM probes	1.5mm	PS101011, PS101012
	2mm	PS100021, PS100022
	3mm	PS100031, PS100032
	4mm	PS100041, PS100042
	5mm	PS100051, PS100052
	7mm	PS100071, PS100072



Medistim PS probe with handle

*Part numbers that end with 1 indicate probes without handle and part numbers that end with 2 indicate probes with handle

Medistim QuickFit™ TTFM probes - PQ probe series

Probe name	Probe sizes	Part numbers*
QuickFit™ TTFM probes	1.5mm	PQ101011, PQ101012
	2mm	PQ100021, PQ100022
	3mm	PQ100031, PQ100032
	4mm	PQ100041, PQ100042
	5mm	PQ100051, PQ100052



Medistim PQ probe with handle

*Part numbers that end with 1 indicate probes without handle and part numbers that end with 2 indicate probes with handle

Medistim Vascular TTFM probes - PV probe series

Probe name	Probe sizes	Part numbers*
Vascular TTFM probe	6mm	PV100061, PV100062
	8mm	PV100081, PV100082
	10mm	PV100101, PV100102
	12mm	PV100121, PV100122
	14mm	PV100141, PV100142
	16mm	PV100161, PV100162



Medistim PV probe with handle

*Part numbers that end with 1 indicate probes without handle and part numbers that end with 2 indicate probes with handle

Medistim Vascular TTFM probes - PA & PB probe series (Note: PA & PB probe series will be phased out by end 2015 and will be replaced with the PV probe series.)

Probe name	Probe sizes	Part numbers*
Vascular TTFM probe	6mm	PB100062
	8mm	PA100081, PA100082
	10mm	PA100101, PA100102
	12mm	PA100121, PA100122
	16mm	PA100161, PA100162



Medistim PA probe without handle

*Part numbers that end with 1 indicate probes without handle and part numbers that end with 2 indicate probes with handle

Medistim Ultrasound Imaging probe - EL probe


Probe name	Probe size	Part number
Medistim L15 Ultrasound Imaging probe	One size	EL100015



Medistim Ultrasound Imaging probe

Medistim Doppler probe - PD probe


Probe name	Probe size	Part number
Medistim Doppler probe <i>Doppler probes are delivered with a removable handle and four stabilizer cups</i>	One size	PD110752



Medistim PD probe

Medistim Cardiac Output TTFM probes - PR probe series

Probe name	Probe size	Part number
Cardiac Output TTFM probe	25mm	PR100251
	20mm	PR100301
	35mm	PR100351



Medistim Cardiac Output probe

APPENDIX B. PROBE CLEANING, DISINFECTION AND STERILIZATION

All Medistim probes are validated for different sterilization methods and should be handled differently according to these procedures.

B.1 General

Cleaning

The surface of the system can be cleaned with hospital grade disinfectants. Avoid excessive use of water or liquids on the system.

For cleaning of probes prior to sterilization, consult the Handling, Cleaning, and Sterilization manual included with each probe.

Disinfection

Approved disinfectants are specified in the Handling, Cleaning, and Sterilization manual included with each probe.

Note that any approved disinfectants are recommended because of their chemical compatibility with product materials, not their biological effectiveness. For the biological effectiveness of a disinfectant, refer to the guidelines and recommendations of the disinfectant manufacturer.

Sterilization

Each probe series is approved for one or more of the following sterilization methods:

- STERRAD® *
- RENO Plasma Sterilizers **
- Ethylene Oxide gas (EtO)
- Steam Autoclave

*STERRAD® is a registered trademark of the company ASP.

** RENO Plasma Sterilizers are manufactured by RENOSEM CO., LTD.

See Appendix A for the sterilization methods approved for each different probe series. Consult also the Handling, Cleaning, and Sterilization manual included with each probe.

It is very important that the corresponding probe cleaning and disinfection instructions are followed carefully.

B.2 Overview of probe cleaning manuals

B.2.1 PA-, PB-, PR-, and PQ-probe series

See Probe Manual for PA-, PB-, PR-, and PQ-probe series PA499001. Handling, Cleaning and Sterilization of the Medistim Probes: Vascular TTFM Probes (PA & PB), QuickFit TTFM Probes (PQ), Cardiac Output Probes (PR).

B.2.2 Doppler probe

See Probe Manual for PD110752 Doppler probe PD499001. Handling, Cleaning and Sterilization of the Medistim Doppler Probe (PD).

B.2.3 PS-probe series

See Probe Manual for the PS-probe series PS499001. Handling, Cleaning and Sterilization of the Medistim QuickFit TTFM Probes, Steam Autoclaveable (PS).

B.2.4 PV-probe series

See Probe Manual for the PV-probe series PV499001. Handling, Cleaning and Sterilization of the Medistim Vascular TTFM Probes, Steam Autoclaveable (PV).

B.2.5 EL15 Ultrasound Imaging probe

See Probe Manual for EL100015 Imaging probe EL499002. Handling, Cleaning and Sterilization of the EL100015 Imaging Probe (EL)



medistim@medistim.com
www.medistim.com

Medistim ASA (Head office)
Økernveien 94
0579 Oslo
Norway
Phone +47 23 05 96 60

Medistim Norge AS
Økernveien 94
0579 Oslo
Norway
Phone +47 23 03 52 50

Medistim Danmark ApS
Gøgetoftens 13
2950 Vedbæk
Denmark
Phone +45 2276 5669

Medistim USA Inc.
14000 25th Ave N. Ste. 108
Plymouth, MN 55447
USA
Phone +1 763 208 9852

Medistim Deutschland GmbH
Bahnhofstr. 32
82041 Deisenhofen
Germany
Phone +49 (0) 89 62 81 90 33

Medistim UK Limited
34 Nottingham South Ind Est
Ruddington Lane
Wilford
NOTTINGHAM
NG11 7EP
Phone +44 (0) 115 981 0871